

```
1  ACLS SQL8.9.x bug-fixes and new features:
2
3  1.  moxyTap-space integration
4      - Completed instuitte staff card, tap to sign in feature to the workspace with
      moxytap-space
5  2.  User profile manager
6      - Account tab: added a quick search tool to pick up the right accounts immediately
      out of large account records
7  3.  Resource manager
8      - Resource QR code: added help tips to editor
9      - Booking alert: added help tips to editor
10 4.  Consumable report
11     - Added full order report by year
12     - Added closed order report by year
13 5.  Message alerts
14     - Included the request type to the message alerts
15 6.  Account manager
16     - Renamed note to description for terminology change
17 7.  Support tab:
18     - Renamed as Requests tab
19     - Renamed Training request to Request for training
20     - Renamed Support request to Request for support
21     - Renamed Sample submission request to Request for sample submission
22 8.  Scan active users who do not set up the default invoiced supervisor on weekly basis
23     - Lab staff receives the weekly alerts
24 9.  Report manager
25     - Sample report: added sample method report and sample type report
26 10. Terms and conditions
27     - Modified terms and conditions PDF access on Dashboard so all institute/lab users
      can access the T&C anytime
28 11. Bugs
29     - Data logbook manager: fixed the retrospective group booking calendar issue
30     - User profile manager: fixed the certificate button to access the editor
      certificate tab issue
31 12. RIMS
32     - Fixed the bugs and added the newly requested features as per modification proposal
33 13. NIF
34     - Fixed the bugs as per modification proposal
35 14. AIIM
36     - Completed the modification as per modification proposal including the cost report
      changes, tracker poster and background image per resource feature
37 15. APPN
38     - Completed the modification as per modification proposal including the special
      weekday and weekend cost model
39 16. Security patch
40     - PHP engine upgrade: 8.4.7
41     - Set PHP expose_php = Off
42
```