

1 ACLS SQL 9.2.x - Bug Fixes & New Features
2
3 1. Supervisor Sign-In
4 + Enhanced supervisor sign-in functionality to provide direct access to their
supervised researcher group, including:
5 - Researcher group profiles
6 - Booking reports for the group
7 - Invoices for the group
8 2. Booking Dashboard
9 - Improved tracker and Tap & Go status user interface for clearer visibility.
10 3. Tracker
11 - Optimised live status tracking scripts for improved accuracy and reliability.
12 4. Calendars
13 + Added filters to properly handle extended ASCII codes in:
14 - Customised booking fields
15 - Booking notes fields
16 5. User Registration
17 - Improved auto-approval functionality to reduce administrator workload.
18 6. MoxyTap Integration
19 - Added automatic sign-out timer feature.
20 - Enhanced MoxyTap Manager user interface.
21 7. Mobile App
22 - Fixed booking permission issue affecting lab spaces enabled for groups of
resources.
23 8. AIIM Customisations
24 + Implemented customisations as per proposal, including:
25 - Operation calendar
26 - Commercial calendar
27 - Service calendar
28 - AIIM finance report
29 9. NIF Customisations
30 - Booking reports: added participant study id
31 10. Database Compliance
32 - Completed compliance upgrades for PostgreSQL v16.9.
33